

## HCR v3 (Choi et al., 2018)

Pre-gaming: Fix samples in fresh 4% PFA in PBS overnight at 4C. Next day, rinse a few times in PBS, and directly embed them in Tissue Tek OCT or NEG-50 and freeze blocks to be stored at -80C (can last for months). Cryosection at desired thickness (usually 14-20 um for zebrafish samples) onto superfrost plus slides and allow slides to air dry in slide box at RT for a few hours to overnight before storage at -80C.

## Day 1:

- 1) Air dry slides for at least 30 minutes
- 2) Fix tissue in fresh 4% PFA in PBS for 10 to 15 minutes
- 3) ProtK in PBS (1 ug/ml: 1 ul into 20 ml) for 5 to 10 minutes depending on tissue thickness. I do 5 minutes for 14 um thick sections
- 4) Rinse in PBS 2x
- 5) Refix in PFA for 5 minutes at RT
- 6) Ethanol dehydration 5 minutes each: 50%, 70%, 100%, 100% Ethanol
- 7) Let slides air dry for 5 minutes at RT
- 8) Put slides in probe hybridization buffer at 37C for at least 10 minutes
- 9) Prepare probe solution (0.4 ul of 2 uM stock HCR probe into 100 ul of probe buffer) (can get away with half the amount of probe)
- 10) Add probe solution to slides and cover in parafilm. Incubate slides at 37C O/N in a moist chamber

## Day 2:

- 1) Remove excess probe by washing 4 x 5 minutes with wash buffer at 37C
- 2) Wash samples 2 x 5 minutes with 5X SSCT (0.1% Tween)
- 3) Put samples in amplification buffer for 30 minutes at RT
- 4) Prepare hairpins separately - 3 ul of 3 uM stock
  - heat at 95C for 90 seconds and cool to RT in a dark drawer for 30 minutes
- 5) Add snap cooled hairpins to 300 ul of RT amplification buffer
- 6) Remove amplification buffer and add hairpin solution. Cover with parafilm and leave in the dark at RT overnight

## Day 3:

- 1) Wash slides in 2 x 5 min in 5x SSCT at RT
- 2) Wash slides 2 x 5 min in 5x SSC at RT
- 3) Dry slide and add antifade mounting reagent and cover slip